



A-44 REGULADOR PARA BAJA PRESION

Regulador para canalizaciones de baja presión. Presión regulable hasta 8 bar. Gases medicinales.

A-44 LOW PRESSURE REGULATOR

Pressure reducer for low pressure pipeline systems. Adjustable outlet from 0-8 bar. Medical gases.



CARACTERISTICAS TECNICAS/TECHNICAL CHARACTERISTICS:

CLASE/CLASS:	IIA
NORMATIVA/STANDARD:	UNE-EN ISO 10524-4
GAS:	GASES MEDICINALES/MEDICAL GASES
CONEXIÓN DE ENTRADA/INLET CONNECTION:	AFNOR/SEO/DIN/FRO/CM/BS
CONEXIÓN DE SALIDA/OUTLET CONNECTION:	1/4M - 9/16M
PRESIÓN DE ENTRADA/ INLET PRESSURE:	≤ 14 bar
ESCALA DEL MANOMETRO/PRESSURE GAUGE SCALE:	0-3 BAR/0-7 BAR/0-16 BAR
CAUDAL ESTANDAR/STANDARD FLOW:	8 Nm ³ /h
DIMENSIONES/ DIMENSIONS:	130 X 100 X 50 mm
PESO/ WEIGHT:	850 gr
CONDICIONES AMBIENTALES DE UTILIZACIÓN/ ENVIRONMENTAL CONDITIONS:	0°C +40°C
CONDICIONES AMBIENTALES DE ALMACENAMIENTO/ENVIRONMENTAL STORAGE CONDITIONS:	-20°C +60°C